

SUSPIRE DISPOSABLE RESPIRATOR **FFP2**

USER INSTRUCTIONS



MAKE SURE YOU HANDS ARE CLEAN AND SANITIZED BEFORE HANDLING THE MASK

- 1) Open the mask with the nose clip at the top and flat it, then pull the ear loops.
- 2) Hold the mask under your chin and then lift to cover your mouth and nose.
- 3) Pull the rings behind the ears with both hands and adjust for the best comfort.
- 4) Using both hands, adjust the shape of the nose pad so that it fits the shape of your nose for the best fit (do not do this with one hand! This will affect the seal of the mask).
- 5) Press the mask onto your face and exhale vigorously. If air escapes from anywhere on the perimeter of the mask, readjust until proper sealing is achieved.

NOTE: The mask can be used for one working shift or multiple times in different days up to 8 hrs in total; in the latter make sure it is stored in a sealed bag and donning and undressing are always performed with clean and sanitized hands.

WARNINGS

- NOT for medical use
- NOT for kids below 10 years old
- NOT for usage in environments where oxygen concentration is below 17 %
- NOT for use in toxic gas environment
- NOT to use with beards or any other facial hairs that interfere with direct touch between the mask and your face skin; carefully shave your face!